Direct Thrombin Inhibitor Assay - 0020302800

Intended use

The HemosIL Direct Thrombin Inhibitor (DTI) Assay is an *in vitro* diagnostic assay for the quantitative determination of the Direct Thrombin Inhibitor Dabigatran in citrated human plasma on the ACL TOP Family of analyzers. The product is intended to be used for the measurement of active Dabigatran levels in situations where measurement is warranted.^{1,2}

Summary and Principle

Direct thrombin inhibitors (DTIs) are anticoagulants that bind to thrombin without the mediation of Antithrombin. Dabigatran etexilate is a prodrug that, once absorbed in the GI tract, is converted to its active form, Dabigatran. Dabigatran is an oral DTI that has been approved for use as an alternative to Warfarin for various therapeutic indications.

The HemosIL DTI assay is a modified (dilute) thrombin time test performed on citrated patient plasma. Citrated patient plasma is diluted in pooled normal plasma (DTI Plasma Diluent - supplied as part of the assay) to reduce interferences from pre-analytical variables. A fixed concentration of reconstituted bovine thrombin is then added to the diluted patient sample, activating the coagulation cascade and converting fibrinogen into fibrin. The presence of Dabigatran in patient samples will have an inhibitory effect on the procoagulant activity of the exogenous thrombin added to the patient sample. The associated clotting time is measured on the ACL TOP Family. A Dabigatran reference curve is plotted on the ACL TOP Family from the clotting time results of known reference plasma standard (HemosIL Dabigatran Calibrators). The concentration of Dabigatran in the patient plasma samples is determined by comparing clotting time values to the reference curve.

Composition

Each Direct Thrombin Inhibitor Assay consists of:

[R]	DTI Thrombin Reagent (Cat. No. 0020302810): 3 x 2 mL vials of a lyophilized bovine thrombin (~1.5 UNIH/mL) with bovine albumin and buffer.
[TD]	DTI Thrombin Diluent (Cat. No. 0020302820): 3 x 2.2 mL vials of a solution containing buffer and heparin neutralizer.
[PD]	DTI Plasma Diluent (Cat. No. 0020302830): 3 x 2 mL vials of a lyophilized human plasma containing buffer, stabilizers and preservatives.

PRECAUTIONS AND WARNINGS:

The Direct Thrombin Inhibitor Reagent contains bovine material. All donor animals were sourced from BSE-free herds. The cattle received ante- and post mortem health inspection by a veterinarian, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious.

The Plasma Diluent uses human material. The human material used in this product was tested by FDA approved test methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious.

The Direct Thrombin Inhibitor Reagent and Plasma Diluent are not classified as hazardous.

Hazard class: None

Hazard statements: None

Precautionary statements: None

Supplemental hazard information:

DTI Thrombin Reagent: Up to 5% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

DTI Thrombin Diluent ≈ 0.5% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

DTI Plasma Diluent: ≈ 100% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

This product is for *in vitro* Diagnostic Use.

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Preparation

DTI Thrombin Reagent – Pipette 2 mL of DTI Thrombin Diluent into a vial of DTI Thrombin Reagent. Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use. Do not shake.

DTI Thrombin Diluent – Ready to use with DTI Thrombin Reagent.

DTI Plasma Diluent – Dissolve the contents of each vial with 2 mL of CLSI Type CLR water or equivalent.³ Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the plasma diluent at 15-25°C for 30 minutes and invert to mix before use. Do not shake.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application/Reference Manual for the complete assay procedure instructions.

Reagent storage and stability

Unopened reagent and diluents are stable until the expiration date shown on the vial when stored at 2-8°C.

DTI Thrombin Reagent – Reconstituted reagent is stable: 48 hours at 2-8°C in closed original vial, or 6 hours on-board the ACL TOP Family at 15°C in the original vial, or 30 days at -20°C in the closed original vial (single thaw only).

DTI Plasma Diluent – Reconstituted plasma is stable: 48 hours at 2-8°C in the closed original vial, or 6 hours on-board the ACL TOP Family at 15°C in the original vial, or 30 days at -20°C in the closed original vial (single thaw only). For optimal stability remove reagents from the system and store them, closed, at 2-8°C, in their original vials.

NOTE: Low relative humidity is associated with increased evaporation of uncapped reagents, which may decrease on-board stability. For optimum on-board stability, laboratory temperature and humidity should be controlled.

Specimen collection, preparation and handling

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to CLSI Document H21-A5 for instructions on specimen collection, handling and storage.⁴

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately:

	Cat No.
Dabigatran Calibrators	0020013400
Dabigatran Controls	0020013500
Cleaning Agent	0009832700
Factor Diluent	0009757600

Quality control

Normal and abnormal controls are recommended for a complete quality control program.⁵ HemosIL Dabigatran Low and High Controls are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8-hour shift in accordance to good laboratory practice. Refer to Westgard *et al* for identification and resolution of out-of-control situations.⁶

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ENGLISH - Insert revision 11/2016

Results

DTI assay results are reported in ng/mL. Refer to the instrument's Operator's Manual for additional information.

Limitations/Interfering Substances

Samples with excessive hemolysis, icterus, or lipemia, should not be used; DTI assay results on the ACL TOP Family are not affected by Hemoglobin up to 300 mg/dL, Triglycerides up to 873 mg/dL, Bilirubin up to 40 mg/dL and Unfractionated / Low Molecular Weight Heparin up to 2.2 IU/mL.

Expected Values

The predictable pharmacokinetic profile of Dabigatran supports a fixed-dose regimen without the need for routine coagulation monitoring. The peak and trough Dabigatran plasma concentration levels have been measured in some clinical trials by the Dabigatran manufacturer and depend on the fixed dosage regime.²

Performance characteristics

Precision

Precision was assessed utilizing 3 lots of plasma on representative members of the ACL TOP Family (ACL TOP 700, ACL TOP 500 CTS, and ACL TOP 300 CTS). Precision was evaluated in accordance with CLSI EP05-A2⁷, for 20 days, with 2 runs per day and 2 replicates per run for each sample level (n=80/instrument/lot), using a specific lot of DTI reagent and both normal and abnormal samples. Test data from a representative instrument, reagent and plasma lot is included below.

ACL TOP Family	Mean (ng/mL)	CV% (Within run)	CV% (Total)
Dabigatran Low Control	45	2.2	3.6
Dabigatran High Control	196	1.8	4.1
Clinical Sample	123	2.0	3.4
Internal Control 1	77	2.8	4.2
Internal Control 2	505	3.5	5.4
Internal Control 3	957	2.3	4.2
Internal Control 4	2035	3.1	4.3

Correlation System:	n	slope (95%CI)	intercept (95% CI)	r	Reference Method
ACL TOP Family	147	1.002 (0.987 to 1.016)	12.01 (5.64 to 18.38)	0.996	Hemoclot Thrombin Inhibitors on ACL TOP Family

A method comparison between the HemosIL Direct Thrombin Inhibitor Assay and the Hemoclot Thrombin Inhibitors (Hyphen Biomed) kit was performed in accordance with EP09-A3, on 147 samples (102 unadulterated Dabigatran patient samples/45 Dabigatran-spiked samples), on a representative member of the ACL TOP Family.

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